


3M™ Hook and Loop fastening system



Introducing HOOK & LOOP fastening tape for medical device construction and closures, the latest component from 3M Medical OEM products. Available with or without adhesive, this fastening system provides a 'low profile' option for the manufacture of medical devices.

SUGGESTED APPLICATIONS: Ostomy device closures; Surgical gowns and other attire; Surgical drape fluid pouches; Tube and cable holders; Monitoring device attachment; Blood pressure cuffs, stretchers, orthopaedic supports and wheel chairs.

The features of this system mean that it is suitable for a wide range of medical devices that require a re-usable and/or adjustable fastener or closure system - make Hook & Loop from 3M your fastening solution!

- **SOFT AND CONFORMABLE** for a comfortable fit
- **STERILISABLE** with ETO and gamma
- **HIGH PEEL STRENGTH** for security during wear and easy adjustment
- **LOW PROFILE** for concealment and comfort
- **LOW LINTING** means reduced secondary attachment
- **3M ADHESION** to films, non-wovens and other substrates

FLEXIBILITY: Work with jumbo or slit rolls, with your converter or a 3M converter, with 3M adhesives or your own nominated adhesive....the choice is yours.



Hook & Loop from 3M comprises of a white low-profile polyolefin hook fastener and a gamma-stable white nylon loop fastener. It contains no natural rubber latex and is available in jumbo or slit roll formats.

PRODUCT DESCRIPTIONS:

- 1. Non-adhesive** - 3M 7330 is a non-adhesive gamma-stable white nylon hook fastener for use with the 7334 white polyolefin hook fastener.
- 2. With Adhesive** - 3M 7331 is a nylon loop fastener laminated to a high tack, synthetic rubber based adhesive designed for use with a low-profile hook fastener (7335), supplied on a two-side siliconized, paper release liner.

For more information, a technical data sheet, samples, or to speak to a customer services representative please call 01509 613267.



Innovative products for a wide range of applications



General Surgical and Medical Devices

Product and function attributes

3M Health Care has a broad experience in the development and production of surgical and medical devices. These can be found in hospitals and clinics around the world, where they have been tried, tested and trusted for years.

3M offers many tape and film components from which to choose when constructing your device - from professional wound care dressings to device attachments to strapping.

Most of our adhesives and tapes are designed to be 'skin friendly' where the devices are attached to patients. Other adhesives and tapes are intended specifically for the construction of a variety of devices made from plastic films, foams, fabrics and treated fabrics.

This comprehensive selection includes components for:

- | | | |
|--|---|-----------------------------------|
| ■ Wound care dressings | ■ Electrodes | ■ Diagnostics |
| ■ Disposable surgical drapes and gowns | ■ Ostomy and continence devices | ■ Cosmetic and skin care products |
| ■ First aid dressings | ■ Personal care and Health Care consumables | ■ Stick-to-skin solutions |

The Trusted Source

For over 50 years, 3M has been a leader within Health Care. During that time, our advanced technologies have contributed to the creation of more than 10,000 new Health Care products worldwide. Working in close co-operation at every stage – from initial product concept through to component supply or actual product design – 3M OEM can link various 3M technologies to product applications and offers a carefully selected line of medical grade tapes, coatings and adhesive technologies, semi finished and finished components together with selected bulk supplies of 3M branded products to medical device manufacturing, design and supply industries.

Supporting you to achieve market success

When it comes to product commercialisation, our expertise in regulatory documentation is second to none. We undertake extensive product testing and provide clinical support packages to enable you to bring your product to market as efficiently and as cost-effectively as possible. What's more, our confidential approach consistently safeguards the integrity of your product development process.

Delivering Innovation, Quality, Service, Value & Confidentiality

For more information, a technical data sheet, samples, or to speak to a customer services representative please call 01509 613267

Please visit us at: www.3m.co.uk/medicaloem



3M Health Care Limited

3M House, Morley Street,
Loughborough,
Leicestershire, LE11 1EP
Phone: 01509 613305
Fax: 01509 613172

3M Ireland

The Iveagh Building
The Park, Carrickmines
Dublin 18, Ireland
Phone: (01) 2803 555
Fax: (01) 2803 509

3M is a trademark of the 3M Company.
© 3M 2009. All rights reserved.