

**3M Medical OEM**  
Silicone Medical Tapes

Bringing your

Products  
to Life

Trusted, Innovative, Global.

**The Powerful Answer  
to Gentle Skin Adhesion**

3M's silicone medical tape is highly processable, patient friendly and available in a variety of formats; making it the ideal solution for stick to skin and stick to device applications.



# Introducing the Powerful Answer

**Medical silicone adhesives** are offering a new gentle solution to applications such as wound care and medical device attachment.



## ...Processable

The variety of backings, formats and liners offers flexibility for device design.



- 3M's medical silicone adhesive tapes are designed for processability, ideal for OEMs and converters.
- Our proprietary method of bonding the adhesive to the backing makes it die-cuttable and reduces potential contamination of equipment.
- Ease of processing with optimum liner release makes this a cost-effective solution for wound care and device fixation applications.
- 3M's absorbent foams and PU films offer a complete solution for advanced wound care products.

|             | Format        | Side A<br>(adhesion)* | Side B<br>(adhesion)* | Carrier   | Colour      | Standard Liner<br>(release)*** | Premium Liner<br>(release)*** |
|-------------|---------------|-----------------------|-----------------------|-----------|-------------|--------------------------------|-------------------------------|
| <b>2475</b> | Single-coated | Si (1.0)              |                       | TPE**     | Transparent | 0.4                            | 0.1                           |
| <b>2476</b> | Single-coated | Si (1.0)              |                       | Non-Woven | White       | 0.4                            | 0.1                           |
| <b>2477</b> | Double-coated | Si (1.0)              | Acr (9.0)             | TPE       | Transparent | 0.4                            | 0.1                           |

\* Adhesion to Steel N/inch 180° peel) \*\* Thermoplastic Elastomer \*\*\* N/inch

To find out how 3M's silicone medical tapes c



# to Gentle Skin Adhesion


## ...Patient friendly



- 3M's breathable and perforated adhesive means long wear-time and good fluid handling capability.<sup>†</sup>
- It delivers a clean and comfortable wear experience, by minimising residue on the skin.
- It's gentle and repositionable, ideal for chronic wounds and compromised skin, and the low profile adhesive means it's conformable to body contours.
- With proven technology from a global health care manufacturer, 3M's catalyst-free adhesive is 'future-proofed' against potential changes in regulations.

## ...Perforated

3M's medical silicone tape delivers consistent, clean perforations due to it's unique adhesive integrity.



### Typical applications

Advanced wound care, monitoring devices, prosthetics, ostomy, tube fixation, skin protection, body tape, neo-natal or geriatric.





**3M Medical OEM UK and Ireland  
3M Health Care Ltd**

1 Morley Street  
Loughborough  
LE11 1EP  
Phone: 0845 8734075  
Fax: 01509 613061  
Email: [medicaloem.uk@mmm.com](mailto:medicaloem.uk@mmm.com)  
Web: [www.3m.co.uk/medicaloem](http://www.3m.co.uk/medicaloem)

**3M Medical OEM Nordic  
3M a/s**

Hannemanns Allé 53,  
DK - 2300 Copenhagen S.  
Denmark  
Phone: +45-4348-0306  
Mobile: +45-2921-6418  
Fax: +45-4396-8596  
Email: [bandersen@mmm.com](mailto:bandersen@mmm.com)  
Web: [www.3m.co.uk/medicaloem](http://www.3m.co.uk/medicaloem)

**3M Medical OEM Germany,  
Austria and Switzerland  
3M Medica**

Hammfelddamm 11  
41453 Neuss  
Phone: +49-2131-14-4457  
Fax: +49-2131-14-4479  
Email: [medicaloem.ge@mmm.com](mailto:medicaloem.ge@mmm.com)  
Web: [www.3mmedica.de/oem](http://www.3mmedica.de/oem)

**3M Medical OEM Belgium  
3M Belgium**

Hermeslaan 7  
1831 Diegem  
Phone: +32-2722-5224  
Fax: +32-2722-5022  
Email: [cstevens@mmm.com](mailto:cstevens@mmm.com)

**3M Medical OEM Netherlands  
3M Nederland B.V.**

Industrieweg 24  
2382 NW  
Zoeterwoude  
Phone: +31-71-5450-450  
Fax: +31-71-5450-688  
Email: [awarmink@mmm.com](mailto:awarmink@mmm.com)

**3M Medical OEM France  
3M France**

Boulevard de l'Oise  
95006 Cergy Pontoise Cedex  
Phone: +33-1-3031-8524  
Fax: +33-1-3031-8376  
Email: [jmpaillaugue@mmm.com](mailto:jmpaillaugue@mmm.com)

**3M Medical OEM Italy  
3M Italy**

Via N. Bobbio, 21  
20096 Pioltello (MI) Italy  
Phone: +39-335-7769-685  
Fax: +39-02-7005-7441  
Email: [lsignorelli2@mmm.com](mailto:lsignorelli2@mmm.com)

**3M Medical OEM Spain  
and Portugal  
3M Health Care**

3M España  
S.A., Juan Ignacio Luca de Tena,  
19 - 25  
28027 - Madrid  
Phone: +34-91-321-6465  
Fax: +34-91-321-6594  
Email: [gdiez@mmm.com](mailto:gdiez@mmm.com)

Visit our web site:  
**[www.3M.co.uk/medicaloem](http://www.3M.co.uk/medicaloem)**